

Early Formulation Screening service

# Accelerating LNP-mediated Nucleic Acid Medicine via a Powerful Ionizable Lipid Library Toolbox

**Dr. Eleni Samaridou**

Early Formulation Screening Service  
Merck KGaA, Darmstadt, Germany


**Millipore**<sup>®</sup>

Expert Pharm/BioPharm  
Products & CTDMO Services

**MERCK**

**Millipore**<sup>®</sup>


Expert Pharm/BioPharm  
Products & CTDMO Services



The life science business of  
Merck KGaA, Darmstadt,  
Germany operates as  
MilliporeSigma  
in the U.S. and Canada.

**MERCK**

Who is  
Merck KGaA,  
Darmstadt, Germany



We are

a vibrant science and  
technology company

Founded in 1668, Merck KGaA, Darmstadt, Germany is the world's oldest pharmaceutical and chemical company and comprises of three unique businesses focused on healthcare, life science and electronics.

## GLOBALLY WE ARE...




**350+**  
HISTORY



**65**  
COUNTRIES

# Our integrated CTDMO offer: Customizable manufacturing from mRNA through lipids to fill & finish

What we do



How we do it


- 1** Quality through Technical differentiation
- 2** Single contact point for a simplified customer experience
- 3** Program de-risk through Regulatory expertise

# The **Early Formulation Screening service** brings in the technical expertise supporting the **preclinical formulation development**



## Preclinical

- **Formulation Screening**
- Formulation Optimization
- Early Process Development
- Analytical characterization & potency testing
- Stabilization and stability testing
- Lipid metabolism & clearance



## Clinical

- Technology Transfer
- Process Scale-up
- Clinical GMP Manufacturing
- CMC document support for regulatory filing




## Commercial

- Commercial GMP Manufacturing
- Aseptic Fill
- Quality Control
- Packaging & Inspection
- Supply Chain

**Early Formulation Screening Service**

**Our integrated LNP CDMO offering**




## Access to a rationally designed ionizable lipid toolbox



### OUR Business model

**IP rights** to top performing ionizable lipids – **cherry-picked from multiple patent families**

**Unique combination**  
of state-of-the-art  
**ionizable lipid libraries**



- ✓ **Preclinical** formulation development (incl. GLP-tox) **under existing Merck licenses**  **no license fees for clients**
- ✓ License required only when ready to enter **GMP phase**

strong lipid candidates identified  
*including performance confirmation in client projects*

- We have **evaluated** internally & together with our lipid partners **hundreds of ionizable lipids** and have **identified the most potent-ones** for your application

Examples of Merck criteria for in-licensing selected ionizable lipids:

- ✓ **Lipid chemistry**, purity, scalability, biodegradability & stability
- ✓ **FTO** through licensing for broad application field
- ✓ **In vivo efficacy & tolerability** – benchmarked to market LNP
- ✓ **Delivery to liver and beyond**

## Screening criteria selected based on customer project

- ✓ Diverse screening criteria can be applied depending on customer project
- ✓ Strong lipid candidates already identified for diverse applications - allowing focused screening panels per application

Example of heatmap generated based on requested formulation screening

Selected LNP	Size	PDI	EE	pKa	Protein Expression				Cell Viability				Hemolysis	Membrane Fusion	Hit	
					HepG2	C2C12	A549	RAW	HepG2	C2C12	A549	RAW				
1																✓
2																✓
3																✗
4																✗
5																✗
6																✗
7																✓
8																✗
9																✓
10																✓
...																
<b>Acceptance criteria</b>	<100 nm	<0.3	>80%	6-7	>BM				≥BM				≤BM	≥BM		✓


PDI: polydispersity index  
EE: encapsulation efficiency

Performance ≥ marketed LNP BM   
Performance < marketed LNP BM 

# Our ionizable lipid toolbox

## Selected lipids from multiple proprietary libraries – benchmarked against marketed LNP

### In vivo tissue tropism tuned by varying the ionizable lipid





**In vivo efficacy ≥ marketed LNP Benchmarks (BM)**



**Ionizable Lipid 1**

**Ionizable Lipid 2**

**In vivo Luciferase activity at 6h post dose (intravenous injection)**



**Ex vivo Luciferase activity at 6h post dose (intravenous injection)**






Legend: Liver (purple), Spleen (pink)



## Prioritizing safety upon administration

### In vivo tolerability study


Hepatic enzymes & cytokine levels measured in blood samples at 6 and 24 hours after injection



**In vivo tolerability comparable to marketed LNP benchmarks**

Our ionizable lipid toolbox


## Lead lipids identified for mRNA vaccine applications – benchmarked against marketed LNP




- ✓ *In vivo* efficacy  $\geq$  marketed LNP Benchmarks (BM)
- ✓ Muscle-to-Liver ratio  $> 1$

# Our ionizable lipid toolbox Strong candidates for mRNA vaccines


## Immunogenicity studies – SARS Cov-2 vaccine




Serum IgG titer



## Challenge Study – SARS Cov-2 vaccine





Percentage of survival of challenged mice



# Our ionizable lipid toolbox

## Validated for multiple nucleic acid payloads

Self-amplifying RNA  
(saRNA)  
delivery




**In vivo saRNA expression**  
 $\geq$  marketed LNP  
**Benchmarks** with  
 tunable pharmacokinetics  
 and peak expression

## Increased efficacy by customizing the formulation composition


**Technical expertise & know-how** to modify the formulation composition in order to **tune the formulation properties, biological performance & tissue tropism**

***In vitro* efficacy**  
increased by ~40times  
by composition optimization



***In vivo* tissue tropism**  
*tuned from liver to spleen*  
by composition optimization\*

**Ex vivo Luciferase activity**  
at 6h post dose (*intravenous injection*)



\*Focusing on the mol% of one specific lipidic excipient

# Formulation Scale-up

## Early-stage process development & in-house technology transfer to large-scale GMP manufacturing


Processes in place for **transition from small to large-scale mixing approaches**

Strong expertise and know-how for efficient **process development – transferable to large scale GMP manufacturing**

Mixing type	Scale increase	Size (nm)	PDI	EE%
<b>Microfluidics</b> (Toroidal mixer)	1-fold	63 ±3	0.18 ±0.04	93 ±5
<b>Microfluidics</b> (High throughput mode / SHM)		67 ±4	0.20 ±0.04	94 ±4
<b>T-junction mixer</b>	100-fold	80 ±4	0.07 ±0.01	94 ±0
<b>Impingement Jet mixer</b>		75 ±1	0.07 ±0.02	97 ±1

SHM = Staggered herringbone mixer  
 PDI = polydispersity  
 EE = encapsulation efficiency

- ✓ **Wide variety of mixing technologies** for different LNP scales to ensure scalability
- ✓ **Manufacturing equipment & processes aligned with our GMP site – technology transfer templates in place**



Potency in selected cell line (normalized)	
<b>Freshly produced LNP</b> (purified by dialysis)	1.00
<b>Frozen LNP</b> (purified by TFF)	1.08

IJM = Impingement Jet mixer  
 TFF = Tangential Flow Filtration

# Lyophilized mRNA-LNP

## No loss of *in vivo* performance upon lyophilization



**Strict ultra-cold chain supply is a significant challenge for RNA-based medicines**  
□ **increased costs & inconvenience for the end-users and for distribution worldwide**



**Lyophilization** of RNA-loaded LNP enables the storage and shipment of LNP formulations at logistically convenient temperatures.

□ *Developing such process remains challenging because of the physical stresses imposed by freezing and drying of the product*



High quality **Emprove®** Excipients for optimal lyoprotectant matrix





Optimized lyophilization **process in place**  
□ *transferable to GMP*



**State-of-the-art** lyophilization & analytical equipment

### **In vivo** performance of lyophilized mRNA-LNP maintained based on our optimized lyophilization process



Confirmed also with siRNA loaded LNP – Data not shown

# Shedding light on metabolism & clearance of novel LNP excipients

New LAUNCH



- ✓ Assess LNP formulations' pharmacokinetics and biodegradability after administration
- ✓ Lipid metabolism and clearance key information related to lipid safety and tolerability
- ✓ Fulfill regulatory requirements to advance to clinical testing

Crucial role of **lipid metabolism/clearance** on the **safety** and **tolerability** of LNP therapies & vaccines

**Regulatory guidance documents** highlight the requirement of **lipidomics data** for novel lipids

moderna  
 > Mol Ther. 2017 Jun 7;25(6):1316-1327. doi: 10.1016/j.jymthe.2017.03.035. Epub 2017 Apr 27.  
**Preclinical and Clinical Demonstration of Immunogenicity by mRNA Vaccines against H10N8 and H7N9 Influenza Viruses**  
 Kapil Bahi<sup>1</sup>, Joe J Senn<sup>2</sup>, Olga Yuzhakov<sup>1</sup>, Alex Bulychev<sup>2</sup>, Luis A Brito<sup>2</sup>, Kimberly J Hassett<sup>1</sup>, Michael E Laska<sup>2</sup>, Mike Smith<sup>2</sup>, Örn Almarsson<sup>2</sup>, James Thomson<sup>2</sup>, Amílcar Mick Ribeiro<sup>1</sup>,  
 Mol Ther Nucleic Acids. 2019 Apr 15; 1b: 1-11  
 Published online 2019 Feb 7. doi: 10.1016/j.omtn.2019.01.013  
**Optimization of Lipid Nanoparticles for Intramuscular Administration of mRNA Vaccines**  
 Kimberly J Hassett<sup>1</sup>, Kerry E Benenato<sup>1</sup>, Eric Jacquinet<sup>1</sup>, Aisha Lee<sup>1</sup>, Angela Woods<sup>1</sup>, Olga Yuzhakov<sup>1</sup>, Sunny Himmelfarb<sup>1</sup>,  
 J Biol Chem. 2018 Oct 10; 293(41):15503-15511. doi: 10.1074/jbc.M118.011010  
**A Novel Amino Lipid Series for mRNA Delivery: Improved Endosomal Escape and Sustained Pharmacology and Safety in Non-human Primates**  
 Staci Sabnis<sup>1</sup>, E. Sathyaith Kumarasinghe<sup>1</sup>, Timothy Salerno<sup>1</sup>, Cosmin Mihai<sup>1</sup>, Tatiana Ketova<sup>1</sup>, Joseph J Senn<sup>1</sup>, Andy Lynn<sup>1</sup>, Alex Bulychev<sup>1</sup>, Iain McFadyen<sup>1</sup>, Joyce Chan<sup>1</sup>, Örn Almarsson<sup>1</sup>, Matthew G Stanton<sup>1,2</sup>, Kerry E Benenato<sup>1</sup>, A. B.  
 Show more  
 + Add to Mendeley Share Cite  
 https://doi.org/10.1016/j.jymthe.2018.03.010 Get rights and content



Molecular Therapy  
 Volume 21, Issue 8, August 2013, Pages 1570-1578  
 Original Article  
**Biodegradable Lipids Enabling Rapidly Eliminated Lipid Nanoparticles for Systemic Delivery of RNAi Therapeutics**  
 Martin A Mainz<sup>1</sup>, Mathiasamy Jayaraman<sup>1</sup>, Shiguo Matsuda<sup>1</sup>, Ju Liu<sup>1</sup>, Scott Barnes<sup>1</sup>, William Querbes<sup>1</sup>, Ying R Tian<sup>2</sup>, Steven M Ansell<sup>2</sup>, Varun Kumar<sup>1</sup>, Jane Qin<sup>1</sup>, Xuemei Zhang<sup>1</sup>, Qianqian Wang<sup>1</sup>, Suni Panesar<sup>1</sup>, Benito Hultobari<sup>1</sup>, Mery Corredo<sup>1</sup>, Julio Hettler<sup>1</sup>, Pochimuthu Kondosomy<sup>1</sup>, Davis Butler<sup>1</sup>, Kalpanasoththi S Rajeev<sup>1</sup>, Bo Pang<sup>1</sup>, Akin Akinc<sup>1</sup>,  
 Anylam  
 PHARMACEUTICALS

- **FDA Liposome Drug Products: Chemistry, Manufacturing, and Controls; Human Pharmacokinetics and Bioavailability; and Labeling Documentation (April 2018)**
- **FDA Guidance for Industry: Safety Testing of Drug Metabolites (Feb 2008)**
- **FDA Draft guidance for Industry: In vitro Metabolism and Transporter-Mediated Drug-Drug Interaction studies (Oct 2017)**
- **EMA: Guideline on the Investigation of Drug Interactions (Jan 2013)**
- **EMA: Reflection paper on the data requirements for intravenous liposomal products developed with reference to an innovator liposomal product (Feb 2013)**

\*not exhaustive list



## Shedding light on metabolism & clearance of novel lipidic excipients in LNP



**Metabolism simulation - in vitro and in vivo metabolism & clearance**

- ✓ Strong expertise and **state-of-art equipment** for accurate quantitation of both **parent lipids & metabolites** (in vitro/in vivo)
- ✓ **Semi-targeted** (via LC-HRMS) & **targeted approaches** (via QqQ MS)

\* LC-HRMS: Liquid chromatography-high resolution mass spectrometry  
 QqQ MS: triple quadrupole mass spectrometry

Powering life-changing medicines

## Your integrated partner for mRNA, lipids and LNP CTDMO Services

Millipore®

- **Best-in class technical support** through **30+ years** of expertise in multiple nucleic acid modalities, lipids and LNPs
- Access to **cherry-picked ionizable lipids** from multiple libraries
- **Large panel of analytical capabilities** (including lipidomics)
- Expertise in LNP **stabilization** and **state-of-art lyophilization** capabilities
- **In-house technology transfer to GMP** manufacturing of LNP

Thank you  
for your attention

We look forward to  
supporting you

Eleni Samaridou, PhD

[Eleni.Samaridou@merckgroup.com](mailto:Eleni.Samaridou@merckgroup.com)



*For more information  
come meet our experts at  
booth #68*



**Technical  
differentiation**



**Simplified customer  
experience**



**Regulatory  
expertise**



The vibrant M, Merck, Emprove and Millipore are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.  
© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.