



# Challenges and Solutions for Modulating Drug Release in OSDF

Raxit Mehta & John Paré

**CRS 2022 Annual Meeting & Expo**

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada

***Advanced Delivery Science***






# Outline


- Emerging innovation trends within the pharmaceutical industry
- OSDF challenges and changes in the developmental road map
- Overcome formulation challenges and achieve target drug release
- Case studies
- Path forward

# Drug Development: Current State




# Innovation Trend: Third Party vs In-house Discovery at Major Big Pharma

Development Work



■ Internal Discovery ■ Discovered by Third Party

Revenue Contribution (total: \$69B)




■ Internal Discovery ■ Discovered by Third Party

Do large pharma companies provide drug development innovation? Our analysis says no  
By Emily H. Jung, Alfred Engelberg and Aaron S. Kesselheim Dec. 10, 2019

# Innovation Trend

## Is Small Pharma Driving Big Pharma Innovation?



- 64% of the 59 FDA approved drugs in 2018 originated from emerging bio-pharma
- Small companies are overwhelmingly driving innovation

Small Pharma Driving Big Pharma Innovation  
<https://www.pharmavoice.com/news/2020-01-pharma-innovation/612330/>

# Small/Virtual Companies



Limited core group of employees responsible for strategic management, regulatory strategy and financial control



Outsource all non-core business functions



Electronic data capture and data submission to regulatory authorities

<https://www.ddw-online.com/the-reality-of-virtual-pharmaceutical-companies-1320-201908/>

# Challenges Faced by the Small/Virtual Companies

<https://www.ddw-online.com/the-reality-of-virtual-pharmaceutical-companies-1320-201908/>

## Goals

- Reduce drug development time
- Reduce R&D cost
- Achieve development activities using minimum resources

## Needs

- Expert opinion/guidance at every stage of the drug product development
- Collaboration with suppliers to mitigate risk
- Regulatory and supply chain support
- Support from CRO/CDMO for development activities

✓ Formulation Expertise with a Trusted Partner

# Formulation Partner Value Contribution




**Expertise of  
Coating, Core  
Formulation and  
Process  
Development**



**Scale-Up and  
Risk Reduction**



**Global  
Manufacturing  
and Support**



**Product Life-  
Cycle  
Management**



**Regulatory  
Support and  
Supply Security**

**Colorcon: Trusted Formulation Partner for over 60 years**



# Services for Speed to Market



**Best-in-Class Products  
for Solid Dosage Forms**



**Global Regulatory  
Expertise**



**Access to Experts and  
Equipment for Trials**




**Processing and  
Characterization  
Capabilities**

# Early Phase Development - Proof of Concept

**Develop prototype tablet or capsule dosage forms when API quantity is limited**


Handles small quantities (mg to grams) of API to provide proof of concept solid oral dosage forms




**API potency requires a higher level of protection**

Validated to safely handle up to Band 4 of 5 APIs (as low as OEL of 1.0 microgram / cubic meter)

# Case Studies



# IR and MR Formulations Across a Wide Range of Dose And Solubility



Formulation: Reference Number 15636 Ver

The following is a suggested starting formulation for your solid oral dose. Colorcon to simplify and reduce the time required for your development information provided. This suggested formula is provided for initial guidance only. Further development and optimization incorporating point formulation only. Further development and optimization incorporating processes and raw materials should be performed. Colorcon technical development activities.

This information has been provided using Colorcon's HyperStart® starting formulation is based upon the following information:

Mr. Jason Hansell  
Colorcon-WP

The HyperStart® starting formulation is based upon the following information:

Ingredient	% w/w	mg / unit
Unspecified Active		
STARCH 1500® (Pregelatinized Starch)	20.00	120.00
METHOCEL™ K100M CR (Hydroxypropylmethylcellulose)	15.00	90.00
Microcrystalline Cellulose (MCC), 90 micron	35.00	210.00
Colloidal Silicon Dioxide	28.50	171.00
Magnesium Stearate	0.75	4.50
<b>Total</b>	<b>100.00</b>	<b>600.00</b>

**Manufacturing Process: Direct Compression**


- Combine the Starch 1500, API and colloidal silicon dioxide; coarsely mix and pass through a 40-mesh (42 micron) screen to break-up agglomerates.
- Add all remaining ingredients, except the magnesium stearate, and blend in a suitable mixer for 15 minutes (or until homogeneously mixed).
- Add the magnesium stearate and blend for an additional 3-5 minutes.
- Perform tablet compression.
- It is recommended to apply a suitable filmcoat to impart mechanical strength, assist packaging, enhance appearance & product stability, and improve patient compliance. Opady® film coating systems available from Colorcon offer the best balance of performance, elegance and productivity. A 3% weight gain of film coating is generally recommended.

**Additional Notes:**

- The actual quantity of active may differ depending upon which salt and/or hydrate form is utilized. In such cases the formulation may be adapted by adjusting the quantity of the diluent(s).
- Tablet hardness may be increased by replacing part of the Starch 1500® with additional microcrystalline cellulose.




# Early-stage Immediate Release Formulation Development



- Formulation that can be converted from capsules to tablets
- Improve blend flow of cohesive drug to reduce weight variation for capsule filling and tableting
- Stabilize moisture sensitive APIs
- Promote content uniformity for low dose APIs

# Capsule to Tablet Conversion to Increase Manufacturing Productivity



Suitable stages of the drug product development: Phase II & III

# Formulation for Direct Compression

Ingredients	Low Dose % w/w	High Dose % w/w
BCS Class I, Very Soluble Drug	5.00	20.00
StarTab® (Directly Compressible Starch)	35.75	28.25
Microcrystalline Cellulose (MCC), PH 102	59.00	51.50
Magnesium Stearate	0.25	0.25
Total	100.00	100.00

Target Tablet Weight  
250 mg for 50 mg dose  
200 mg for 10 mg dose

# Capsule to Tablet – Key Challenges

Effect on  
drug release


Polymorph  
change

Flowability







# Drug Release Evaluation




Criteria: >85% drug release  
in 15 min or less

→ No additional in vivo study

# Coating of Tablets: Overcoming Concern of Polymorph Change



# Blend Flowability



Good blend flow for low dose formulation and passable flow for high dose formulation

# Capsule to Tablet – Key Challenges

Effect on  
drug  
release


Polymorph  
change

Flowability


✓ Successfully addressed the challenges



# Low Dose and Directly Compressible Formulations



# Drug Attributes



*Total tablet weight was 400 mg; 9.5 mm round std concave tablet  
Direct compression using rotary tablet press @ 35 rpm; 15 kN compression force*


# Naproxen Tablet Formulation

Ingredients	Formula A (StarTab)	Formula B (Lactose)
Naproxen (micronized)	1.00 %	1.00%
StarTab® (Directly Compressible Starch)	98.50%	-----
Lactose Monohydrate (FastFlow 316)	-----	96.50%
Crospovidone	-----	2.00%
Magnesium stearate	0.50%	0.50%
<b>Total</b>	<b>100.00%</b>	<b>100.00%</b>

*Total tablet weight was 400 mg; 9.5 mm round std concave tablet*

*Direct compression using rotary tablet press @ 35 rpm; 15 kN compression force*

# Dissolution Comparison




*Dissolution: USP monograph (phosphate buffer pH 7.4; 900 mL, Paddle, 50 rpm)*



# Blend and Content Uniformity

Formula	Description	Description	Blend Uniformity, %	Tablet Content Uniformity, %
Formula A	StarTab ; DC	Top right	95	----
		Top left	96	----
		Bottom	94	----
		Composite	94	----
		Tablets, 5 min	----	97.23 ± 0.34 *
		Tablets, 10 min	----	97.32 ± 0.60 *
Formula B	Lactose; DC	Top right	84	----
		Top left	87	----
		Bottom	88	----
		Composite	87	----
		Tablets, 5 min	----	86.50 ± 0.63 *
		Tablets, 10 min	----	86.06 ± 0.86 *

\*Mean ± % RSD



# Low Dose and Directly Compressible Formulations

- Formulation based on StarTab® resulted in improved blend and content uniformity
- StarTab is recommended for low dose formulations to improve content uniformity after simple blending and direct compression



# Bilayer and Fixed Dose Combination (FDC) Technologies for Biphasic (IR+ER) Performance



- Bilayers for fixed dose combinations
- Multiparticulates, including mini-tablets
- Coating expertise for layering of low dose drugs on tablets

# Bilayer Tablet Composition



**IR layer**

Direct Compression		
Ingredients	% w/w	mg/tablet
Glimepiride	0.500	1.000
STARTAB®	33.080	66.160
Lactose monohydrate	65.670	131.340
Dye color Iron Oxide yellow	0.500	1.000
Magnesium Stearate	0.250	0.500
TOTAL	100.000	200.000

Hardness = 22 - 25kP

**ER layer**

Granulation		
Ingredients	% w/w	mg/tablet
Metformin Hydrochloride	71.400	499.800
METHOCEL K100M Premium	23.300	163.100
MCC	4.000	28.000
Colloidal Silicon Dioxide	0.900	6.300
Magnesium Stearate	0.400	2.800
TOTAL	100.000	700.000


Hardness = 18- 22kP

**Tablet dimensions: 16.5 x 8 mm, caplet**  
**IR + ER = 900 mg**


# Bilayer Tablets - Release Profiles

## *Uncoated and Coated tablets*

Release Profiles of (A) IR Glimepiride and (B) ER Metformin HCl from Uncoated and Coated Bilayer Tablets



900 mL pH 7.8 Phosphate buffer + 10% SLS, using paddle with sinkers @100 rpm for 90 mins, analyzed using HPLC @228nm



USP Apparatus II with sinkers @100 rpm,  
1000 mL phosphate buffer pH 6.8, analyzed using UV @233nm

# Bilayer Tablets - Release Profiles

*Uncoated and Coated tablets*



Uncoated and coated tablets showed similar dissolution profiles



> 75% glimepiride released in 10 min (IR) and > 90% metformin was released in 8h (ER)



# Film Coating of Low Dose, Low Soluble API on ER Matrix Core

## Formulation Development Activities

- ER matrix core: metformin hydrochloride, aqueous granulation
- IR film coat: glimepiride, aqueous based coating system

## Quality evaluation parameters for film coating layer

- Content uniformity, Assay and dissolution testing of glimepiride



# Composition of Film Coated Tablet

ER Matrix Core  
Metformin hydrochloride

ER Matrix Tablet	Granulation		
	Ingredients	% w/w	mg/tablet
	Metformin Hydrochloride	62.50	500.00
	METHOCEL K100M Premium	32.20	257.60
	MCC	4.000	32.00
	Colloidal Silicon Dioxide	0.900	7.20
	Magnesium Stearate	0.400	3.20
	TOTAL	100.000	800.000

*18 x 8 mm caplet shape*

IR Film coating  
Glimepiride

Component	Functionality
Glimepiride, $d_{90} = <10\mu\text{m}$	Active
Opadry® II	Binding agent
SLS	Wetting agent
Water	Dispersion medium



# Film Coating Optimization Trials

Components	%		
	Trial 1	Trial 2	Trial 3
Glimepiride	2.50	1.25	0.833
Opadry® II	96.50	97.75	98.167
SLS	1.00		
Total, %	100		
% Weight gain	5	10	15
Final tablet weight, mg	840	880	920

# Results - Glimepiride IR Layer


Glimepiride  
Assay %

Trial 1 5% WG	Trial 2, 10% WG	Trial 3, 15% WG
112%	99%	106%
		


- Coated tablets appeared smooth and uniform
- Use of dilute dispersion for drug layering application helped meet assay specifications of low dose API, glimepiride

# Results - Glimepiride IR Layer

## Content Uniformity



## Release Profile



USP Apparatus II, 100 rpm, 900 mL pH 7.8 Phosphate buffer + 1% SLS


- Dilute drug layering dispersion provided uniformity of API distribution with AV of 9.69
  - impact of higher WG ~ use of dilute API dispersion
- Low dose, low soluble API like glimepiride meets dissolution specs
  - use of wetting agent (SLS) and increased ratio of Opadry: drug helped improve dissolution of glimepiride

# Your Formulation Partner

With different pressures in the development of new pharmaceutical products:



**Speed To  
Market**



**Access to  
Expertise**



**Changing  
Legislation**

...it's important to choose the right collaboration partner

# Colorcon Supports Your Early Stage Formulation Development



Reduce  
Time to Market



Lower  
Total Cost



Make a Positive  
Contribution to  
Quality



Be a  
True Global  
Partner

Colorcon minimizes risk by providing the right solution for your business needs

# Please stop by our booth to learn more:

Booth #217



Let's connect and collaborate to discuss your formulation needs and to achieve desired release profile

[rmehta@colorcon.com](mailto:rmehta@colorcon.com)

[jpare@colorcon.com](mailto:jpare@colorcon.com)