

# Exploring Promising Alternatives to In Vitro Performance Testing of Dry Powder Inhalers

Ken Boda  
Dissolution Product Specialist  
Agilent Technologies

# Overview



Overview of Dissolution



Pulmonary Delivery and Challenges in Release Rate Testing



Case Study I using an Enhancer Cell



Case Study II using Small Volume Apparatus 7 (400-DS)



Promising option using NanoDis – IR and ER Case Studies



Summary & Key Takeaways

# What is Dissolution?

Dissolution is the Determination of a drug's rate of release In Vitro under a specific set of conditions.

This release rate is correlated to in vivo data and can be used to determine the acceptability of a formulation for use.



# The Importance of Dissolution

Of all the tests performed on pharmaceutical products, dissolution is the only one that determines the ***performance*** of a formulation.

Dissolution is essential in making sure that drug is administered in a **safe and effective** manner at all stages.

# What do we test?

Dissolution is not just about orally ingested products. It is used for tablets, capsules, powders, ointments, creams, gels, osmotic pumps, transdermal patches, implants, bone cements, ocular implants, drugs on chips, medicated contact lens, stents, coated beads, catheters, wound care products, etc....



Tablets/capsules




Suspensions




Ointments




Transdermal patches




Pacemaker Lead



Contact Lens



Ophthalmic Implant



Drug Coated Stents

# Dry Powder Inhalers

Inhalation Route of Administration offers many advantages for drug delivery:

- Non-invasive
- Avoids 1<sup>st</sup> pass liver metabolism
- Large epithelial surface area with good permeability
- Direct transport into the bloodstream allows optimal systemic delivery



# Dry Powder Inhalers

- Particle sizes are typically sized 0.5 – 4.0 $\mu$ m
  - <0.5 risks being exhaled
  - >4 can get trapped in the brachial airways and cleared through coughing or swallowing
- Prior to Dissolution Testing, proper fractionation of the dose from the inhaler should be obtained.
- USP <601> Inhalation and Nasal Drug Products: Aerosols, Sprays, and Powder – Performance Quality Test discusses particle evaluation equipment




# Dissolution Attempts with Dry Powder Inhalers

Dissolution Testing with Dialysis Chambers has been attempted for DPIs but there are significant issues with this approach:

- Insufficient mixing within the dialysis chamber itself
- Permeation through the membrane becomes the rate limiting step
  - Sink condition may not be maintained
  - Critical burst phase may not be captured



Float-A-Lyzer with holder and evaporation cover  
(p/n 12-2066).





# Promising Options for Metered Dose Inhaler Dissolution

- Enhancer Cell for use with Modified USP Apparatus 2
- USP Apparatus 7 Small Volume System (400-DS)
- NanoDis nanoparticle filtration system for use with USP Apparatus 1 and 2



# Case Study I: Dissolution of Fine Particle Fraction from Truncated Anderson Cascade Impactor with an Enhancer Cell

In this study, the collection of fine particle fractions (FPF) utilized a special PTFE funnel over the final small collection plate (sCP) contained in the truncated ACI.

After determining the FPF for the drugs under study, various stages were removed to truncate the ACI which allowed a simplified and robust means of collecting only the fine particles on the final sCP.

The sCP was coated with a film of silicon oil to retain the particles. Next, a portion of adhesive tape equivalent in size to the 4.0 cm<sup>2</sup> Enhancer Cell was used to recover the dose deposited on the sCP and mounted in the Enhancer Cell for testing facing upward.




Fig. 1. Schematic diagram of truncated ACI set up with the cross-section of the additional insert on stage F for the collection of FPF on single stage.

# Case Study I: Dissolution of Fine Particle Fraction from Truncated Anderson Cascade Impactor with an Enhancer Cell



- Described in USP <1724> Semisolid Drug Products as Immersion Cell Model A
- Works in a traditional Apparatus 1 and 2 system
- Typically used for testing semisolids
- Good alternative to Franz Cell testing, allows use of existing systems
- Available in 0.5, 2, and 4cm<sup>2</sup> openings
- Membrane selection dependent on application

# Case Study I: Dissolution of Fine Particle Fraction from Truncated Anderson Cascade Impactor with an Enhancer Cell

## Method Parameters:

- 50 RPM Apparatus 2 with mini-paddle
- 200mL mini-vessel set up with Enhancer cell
- Media: Gamble's Solution with 0.2% v/v Tween 80 to mimic simulated lung fluid at 37°C
- Paddle placed 1cm above enhancer cell
- 2mL samples taken with media replacement



# Case Study I: Dissolution of Fine Particle Fraction from Truncated Anderson Cascade Impactor with an Enhancer Cell




Fig. 4. Dissolution profile of FP (×), SX (●) and SS (▲) in mSLF.

# Case Study I: Dissolution of Fine Particle Fraction from Truncated Anderson Cascade Impactor with an Enhancer Cell

## SUMMARY:

The Method was able to differentiate between 3 particle formulations successfully.

- SS is BCS Class 1 with high solubility as well as high surface area
- SX which is highly soluble showed a rapid burst profile
- FP follows a logarithmic function

These were expected given their relative solubilities. In vivo studies support the slower release of glucocorticoids such as FP.

**Use of the Enhancer Cell allowed for a discriminatory dissolution method using commercially available equipment already present in the USP, with a slight modification.**

Use of this approach could allow for a reproducible means of testing these products and would also be adaptable for automation such as autosamplers and/or online UV-Dissolution analysis.

## Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology

- Drug Loaded Pellets were made through extrusion-spheronization
- 25% w/w MMC, 25% w/w lactose, and 50% w/w metformin was used
- Drug pellets were then transferred to a fluid bed coated and sustained release pellets with 8% and 10% w/w coat weight gains were created
- Pellets were then dried at 60 C for 12 hours

322

*M. Xu et al./International Journal of Pharmaceutics 478 (2015) 318–327*

**Table 2**

Physicochemical properties of sustained release pellets with 8% and 10% coats.

Coat (%)	$D_{50}$ ( $\mu\text{m}$ )	$SP_{\text{pel}}$	Aspect ratio	Roundness	$W_s$ ( $\mu\text{g}$ )	Drug content (%)	Coat thickness ( $\mu\text{m}$ )
8	$813 \pm 24$	$0.35 \pm 0.01$	$1.11 \pm 0.05$	$1.13 \pm 0.01$	$516 \pm 10$	$48.7 \pm 4.0$	$13.9 \pm 2.2$ (15.7%*)
10	$849 \pm 56$	$0.40 \pm 0.03$	$1.10 \pm 0.05$	$1.12 \pm 0.01$	$513 \pm 8$	$49.4 \pm 3.5$	$15.7 \pm 0.9$ (5.9%*)

\* The relative standard deviation of pellet coat thickness.

# Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology

Initial Studies were performed with USP Apparatus 1 and 2. Data showed high variability and different phases of dissolution release.

This data is explained by floating and clumping behavior of the beads in both apparatus types.




Fig. 3. Cumulative drug release profiles of ensemble of sustained release pellets with 8% and 10% coats using USP dissolution apparatus 1 (basket) and apparatus 2 (paddle).



# Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology

## USP Apparatus 7 – Small Volume – 400-DS


- Originally designed for medical device market
- Low dose, high potency drugs
- Essentially no evaporative loss
- 3 - 12mL test conditions



# Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology

## Small-volume USP App 7: 400-DS

- Sample holder reciprocates 2 cm per USP <724> Drug Release
- Mixing from Reciprocation avoids the challenges of agglomeration of particles as well as the challenges of floating
- Product is completely submerged and free-flowing within the dissolution media



# Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology

## Holder Types for the 400-DS

Products that can work well with this method include:

- Subcutaneous implants
- Stents
- Contact Lenses
- Wire/Tubing/Leads
- Birth Control (IUDs, Implants, etc.)

Product must be **non-disintegrating** due to no filtration.



# Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology

## 400-DS Apparatus 7 Study

In this study, particles were separated into various fractions using an Anderson Cascade Impactor as with the previous study.

The fractions were collected onto a filter, and another filter was placed on top. This particle/filter sandwich was then rolled tightly and placed within the 400-DS PEEK basket holder



## Case Study II: Evaluation of the coat quality of sustained release pellets by individual pellet dissolution methodology




Fig. 4. Cumulative drug release profiles of individual sustained release pellets with 8% and 10% coats and their linear correlation values.

- 400-DS method parameters were optimized with changes to dip speed, dip interval, and dissolution media volume
- Linear dissolution release rates were able to be determined through the optimized methods with R2 values of 0.991 and 1.000 for the 8% and 10% coated products, respectively
- More variability was seen in the 8% v/v and indicates a greater variability of coat quality compared to the 10%.

# Promising Options Using the NanoDis System

## Cross Flow Filtration for Nanoparticles




# NanoDis - Hollow Fiber Filters

## Cross Flow Filtration for Nanoparticles

- Used for manufacturing for purification and solvent removal
- Tangential flow of the filtration medium “different from the dead-end filtration”
- No filter cake
- No blockage of the filters

Dissolution Media  
with Dispersed API  
loaded Nanoparticles  
+ Dissolved (free)  
API




Dissolution Media  
with Dispersed API  
loaded Nanoparticles

# NanoDis - Hollow Fiber Filters

## Cross Flow Filtration for Nanoparticles

- Nanoparticle formulations have great potential for pharmaceutical development
  - especially for low soluble substances
- Small particle size of nanoparticles results in
  - greater surface area
  - greater bioavailability of the Active Pharmaceutical Ingredient (API)






# NanoDis - Hollow Fiber Filters

## Cross Flow Filtration for Nanoparticles

However, conventional dissolution techniques present significant challenges:

- difficult to filter using conventional membrane technology, thus the dissolution process continues
- filters may rupture and dispersed nanoparticles pass into the measuring vial
- nanoparticles tend to clump within dialysis membranes which can become blocked and incorrect results are observed
- unsuccessful separation of nanoparticles from the dissolution medium
- negative impact on the ability to accurately predict the true API dissolution profile and selection of the best lead formulation



# The NanoDis System – what is it?

A workflow that ensures *in vitro* dissolution of nanoparticle API



# NanoDis System



Unique to Agilent, the NanoDis System addresses the pharmaceutical industries need for compliant, semi-automated testing of nanoparticle drug formulations.

The NanoDis System enables formulation scientists to produce accurate release profiles of API's utilizing conventional dissolution apparatus:

- A market first that combines Agilent's conventional USP App 1 & 2 (708-DS) with automated sampling (850-DS), cross-flow filtration, and flexible software
- An automated, compliant solution ideal for regulated environments


# The NanoDis System – what is it?

A workflow that ensures *in vitro* dissolution of nanoparticle API





# The NanoDis System – what is it?

A workflow that ensures *in vitro* dissolution of nanoparticle API



# The NanoDis System – what is it?

Powered by flexible, compliant Dissolution Workstation Software



Method Editor

Parameters | Dissolution | Sample Timepoints | Prefill | Notifications | Auto Export | NanoDis

Sample Information

Product Name: Prednisone

Notes: For PQ

User Defined Labels

Label 1: Expiration Date

Label 2: Lot #

Label 3: Group

Cannula Filter Type: 70 micron Full Flow

Sampling Parameters

Sample Volume: 7 mL

Prime Volume: 5 mL

Purge Volume: 5 mL

Active Channels: 6

Media Replacement

Replacement Volume: 0

Waste Drop Vol: 0 mL

Media Volume: 900 mL

Media Replacement

Replacement Volume: 0

Waste Drop Vol: 0 mL

Media Volume: 900 mL

Media Replacement

Replacement Volume: 0

Waste Drop Vol: 0 mL

Media Volume: 900 mL

Error Tolerance (°C)

Temperature: 0.5

Speed: 2 %

Profile Interval: 00:00 HH:MM

Clean System After Method Run



Enable: Vol: 10 mL Cycles: 1

Change Management

Restrict method editing rights to the method owner

Owner: Everyone

Cancel OK



Method Editor

Parameters | Dissolution | Sample Timepoints | Prefill | Notifications | Auto Export | NanoDis

Properties

Bath Temperature: 37.4 °C

Vessel Temperature: 37.0 °C

Spindle RPM: 50

Final Spin RPM: 250

Final Spin Duration: 00:00 MM:SS

Media Volume: 900 mL

Apparatus Type

Paddles

Sequential Start: 00:10

Vessel Table Level

Enable Pre-Test Level Verification (requires 280-DS IM)

Tolerance: 0.5

Operation Checks

Enable Pre-Test Verification of Apparatus Components and Environment

Vibration Monitoring



Enable Vibration Monitoring (requires 280-DS Instrument Module)

Frequency Bandwidth of Interest: 0 to 250 Hz

Transient Event Threshold: 30.0 mG

Persistent Event Threshold: 10.0 mG over 05:00 MM:SS

Cancel OK



Method Editor

Parameters | Dissolution | Sample Timepoints | Prefill | Notifications | Auto Export | NanoDis

Enable NanoDis

NanoDis Method Properties

Pre-Test Filter Conditioning

Peristaltic Flow-Through Duration: 20 s

Syringe Purge Volume: 5 mL

Peristaltic Air Purge Duration: 20 s

Timepoint Sampling Properties

Filter Outer Cylinder Rinse Volume: 10 mL

Filter Outer Cylinder Rinse Cycles: 1

Peristaltic Filter purge duration: 45 s

Purge filter toward sample cannula

Pre-Timepoint Filter Conditioning

Peristaltic Flow-Through Duration: 20 s

Syringe Purge Volume: 5 mL

NanoDis Filter Type: Generic Filter

Cancel OK



# The NanoDis System – what is it?

Using the right methodology is critical in drug development and quality control of nanoparticle formulations.

Critical steps:

1. complete separation of API from the nanoparticle
2. a quick separation method - especially for immediate release formulations
3. a blockage & rupture free filtration solution
4. quantification of only the dissolved API



# The NanoDis System – what is it?


- **Separation** of nanoparticles from dissolution media independent of particle size
- **Online** and simultaneous particle separation
- **Accurate** *in vitro* dissolution of nanoparticles and quantification of only dissolved API
- Ability to test immediate and extended-release nanoparticle formulations
- A system that is free from blockages and filter ruptures






# The NanoDis System – what is it?

- An uninterrupted dissolution process as media and any dispersed nanoparticles are recycled back to the 708-DS vessels
- **Rugged and robust** methodology using conventional USP apparatus. Methods are easily transferred to QA/QC.
- **Automatically** maintained and consistent fill-volumes during analysis
- Full **compliance** - all methods are operated, electronically documented and auditable (DWS)
- **Peace of mind** knowing that a fully compliant Dissolution Testing solution manufactured to stringent quality control procedures has been implemented




# NanoDis System > Filter Efficiency Test



*Separation of nanoparticles from dissolution medium after sampling*

# NanoDis System > Filter Efficiency Test




No difference between unfiltered sample, 0.45 µm, 0.2 µm filtered sample

**COMPLETE PERMEATION OF NANOPARTICLES THROUGH FILTER**

Nanoparticles are absent in Nanodis samples

**COMPLETE SEPARATION OF NANOPARTICLES WITH NANODIS**




All dissolution samples are colorless with Nanodis

# Dissolution of Immediate Release (IR) Nanoparticles

# Immediate Release Nanoparticles

## Case Study



Ibuprofen  
nanoparticles



### Selection of dissolution medium

- A dissolution medium in which nanoparticles are 100% dispersible
- A dissolution medium where the API is limited soluble




### pH 4.5 Acetate buffer

- Particles do not agglomerate
- Solubility limited in the medium
- The aim of dissolution is to determine
  - ✓ ONLY the dissolved API, and
  - ✓ NOT dispersed nanoparticles

# Immediate Release nanoparticles

## Case Study



Ibuprofen  
nanoparticles in pH  
4.5 acetate buffer  
200 mg in 900 mL

**Nanoparticles  
dispersed in  
dissolution  
medium (60%)**

**Nanoparticles  
dissolved in  
dissolution  
medium (40%)**

**FAILURE to separate  
nanoparticles from dissolution  
medium...**


- post-sampling analysis with HPLC shows 100% dissolution
- because dispersed samples dissolve during HPLC analyses

**SUCCESSFUL separation of  
nanoparticles from dissolution  
medium...**

- post-sampling analysis with HPLC shows max 40% dissolution

# Immediate Release Nanoparticles

## Case Study



- 100% Dissolution is reached with **0.45 µm** filter although the API solubility is limited to 40%. This shows the complete permeation of nanoparticles through the 0.45 µm filter
- Lower dissolution values were observed for **0.2 µm** filter due to the particle size of Ibuprofen where a significant amount of particles are larger than 220 nm
- **0.02 µm** filter was blocked during sampling resulting in lower and inconsistent dissolution profile

### NanoDis System

- Only dissolved API was characterized using NanoDis enabling true *in vitro* dissolution testing.
- Nanoparticles were completely retained by the filter

# Dissolution of Extended-Release (ER) nanoparticles



# Extended-Release Testing

## Case Study

**Rhodamine PLGA  
Nanoparticles**



Dialysis bag

Nanodis


**Dissolution studies**



Quantification of  
Rhodamine with HPLC

# Extended-Release Testing

## Case Study



### NanoDis System

- quantifies the real dissolution kinetics for the extended-release formulations
- enables true *in vitro* dissolution testing.

Dissolution through dialyses membrane is controlled by

- permeation coefficient, and
- supersaturation in the dialyses membrane

# Summary & Key Takeaways

- Currently a lack of clear guidance for Dissolution Testing of Metered Dose Inhalers
- Many attempts with traditional filtration and dialysis techniques have severe challenges associated with them
- There are new promising approaches to these types of products using:
  - NanoDis Cross-Flow Filtration
  - Enhancer Cell Immediate Cell
  - 400-DS Small Volume Apparatus 7



# Summary & Key Takeaways

## Acceptable Method Requirements



Low-to-moderate variability

Complete Release (85%+ or Asymptote)

Characterizing timepoints below 85%

Challenged with other formulations

Evaluated to ensure they are rugged and reproducible enough for repeated testing by multiple people/units

# Method Development Resources

- USP <1092>
- FDA Dissolution Method Database (<http://www.accessdata.fda.gov/scripts/cder/dissolution/index.cfm>)
- Dissolution Discussion Group ([www.dissolution.com](http://www.dissolution.com))
- Dissolution Technologies ([www.dissolutiontech.com](http://www.dissolutiontech.com))

# Questions?

Contact me at:

- [Ken.Boda@Agilent.com](mailto:Ken.Boda@Agilent.com)
- [Dissolution.Hotline@Agilent.com](mailto:Dissolution.Hotline@Agilent.com)
- [www.dissolution.com](http://www.dissolution.com)
- or find me on LinkedIn

