

Simplifying Progress



LipidBrick®: An innovative lipid library for efficient and tunable mRNA-LNPs

CRS, July 2024

SARTORIUS

Delivery of RNA and DNA is not one size fits all

- Lipid nanoparticle solutions have shifted to a mostly one size fits all format:
 - Nearly identical lipid compositions and biodistributions in vivo
 - Lack of formulations tailored to specific payload and applications
- The issues with a one size fits all approach are abundant:
 - Lack of tissue specificity and/or broad biodistribution
 - Applying complex API to the needs of the excipient
 - Lower efficacies than needed

Theranostics: Yasmin Granot-Matok et, al. pp3397 ISSN1838-7640




Unique Formulations Solve Current Challenges

1	Biodistribution	Broaden the Approach	<ul style="list-style-type: none">▪ Current LNPs accumulate in the liver▪ Need for extrahepatic delivery to broaden therapeutic applications
2	IP	Freedom to Operate	<ul style="list-style-type: none">▪ The current IP landscape for LNPs is complex▪ Components and formulations are both patented▪ Need for secure IP to ensure translation to clinic
3	Efficacy	Uptake and Response	<ul style="list-style-type: none">▪ Need for safe, non/low immunogenic delivery solutions that can be taken up abundantly▪ Need for high transfection of cells
4	API Driven	One Size ≠ Fit All	<ul style="list-style-type: none">▪ Optimization of delivery systems for specific RNA▪ Ability to design formulations that are unique


Focus on: Intellectual Property

- Patents and patent infringement remain extremely relevant topics in 2024
- Majority of preclinical research is done with MC3 (Onpattro), SM-102 (Moderna), ALC-0315 (Acuitas)
- Lack of clear pathway to the clinic for several lipids




<https://www.nature.com/articles/s41587-021-00912-9>

Deep History in Nucleic Acid Delivery




Polyplus has developed multiple GMP products

Simplifying Progress in DNA | RNA Therapeutics



Critical Raw Materials



Services



Bioreactor Technologies




Lipids




Transfection Reagents




pDNA




recHuman Albumin



pDNA Engineering



pDNA Manufacturing



PD Services



Multi-parallel Bioreactors



Scalable Stirred-Tank Bioreactors

End-to-End Solutions Supporting Scalable Nucleic Acid-based Therapeutics



Bio Analysis



Data Analytics



Fluid Management



Downstream Processing



At-Line HPLC



Analytical Monoliths



DOE and MVDA Software



Sterile Connection & Disconnection Systems



Single Use Bags



Storage & Mixers



TFF Systems & Consumables



Filters




Monoliths




SU Batch Chromatography

ATS Solutions Across the pDNA | mRNA Workflow




ATS Solutions Across the pDNA | mRNA Workflow



LipidBrick® library


LipidBrick® - A breakthrough technology

- New generation of cationic lipids to form LNPs
- Based on an **IMidazolium** polar head, a heterocycle, leading to **reduced toxicity** as compared to traditional cationic lipids
- IM21.7c has been **used and validated for years** in different types of formulation:
 - jetMESSENGER® (lipoplexes) → *in vitro* mRNA transfection
 - *in vivo*-jetRNA®+ (liposomes) → *in vivo* & *ex vivo* mRNA transfection
- **Sole proprietary technology of Polyplus®** (Now part of Sartorius)
 - EP18306417 (Granted in EU, China, Hong Kong, Japan. Pending or under review in the US and other countries)
 - WO2024008967A1 (Published in Jan 2024)
 - FTO analysis (US & EU patents) performed by ext. IP attorney




 **Cationic lipid**
= Fixed positive charge

DOTMA - Previous generation cationic lipid



Applying LipidBrick® to Conventional LNP Design



Active substance = Nucleic acids



mRNA


Carry genetic information to be translated, *in fine*, into a protein of interest

Excipient/Carrier = LNPs (lipid nanoparticles) composed of a mixture of lipids



Cationic / ionizable lipid

Active positively charged lipids will interact with negatively charged nucleic acids and cell membrane (improve cell uptake)



 **LipidBrick®**



Helper lipid

Structural lipids provide LNP rigidity and stability



PEGylated lipid


Provide steric stabilization and prolong blood circulation





Cholesterol

Increase membrane fluidity and stability, promote membrane fusion




Design of innovative lipids







Design of innovative lipids








Design of innovative lipids



Design of innovative lipids




Design of innovative lipids



LipidBrick® - A Full Library Offering Plenty of New Options

IM21.7c
CAS n°: 2416939-42-7




<p>IM16c</p> <p>Chemical structure of IM16c, a lipid molecule featuring a long hydrocarbon chain (black) and a branched alkyl chain (black) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>	<p>IM22c</p> <p>Chemical structure of IM22c, a lipid molecule featuring a long hydrocarbon chain (red) and a branched alkyl chain (black) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>	<p>IM25c</p> <p>Chemical structure of IM25c, a lipid molecule featuring a long hydrocarbon chain (red) and a branched alkyl chain (black) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>
<p>IM13c</p> <p>Chemical structure of IM13c, a lipid molecule featuring a long hydrocarbon chain (black) and a branched alkyl chain (blue) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>	<p>IM12c</p> <p>Chemical structure of IM12c, a lipid molecule featuring a long hydrocarbon chain (red) and a branched alkyl chain (black) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>	<p>IM3c</p> <p>Chemical structure of IM3c, a lipid molecule featuring a long hydrocarbon chain (black) and a branched alkyl chain (black) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>
<p>IM15c</p> <p>Chemical structure of IM15c, a lipid molecule featuring a long hydrocarbon chain (black) and a branched alkyl chain (black) connected to a pyrazole ring. The pyrazole ring is substituted with a methyl group (Me) and a butyl group (n-Bu). The molecule is shown as a chloride salt (Cl⁻).</p>		

Proof of Concept: LipidBrick[®]-based LNPs



LipidBrick® are versatile

- LipidBrick® IM21.7c has been shown to be suitable with:
 - Different compositions:
 - As complete substitution of ionizable lipid → 4-lipid LNP
 - In addition to ionizable lipid → 5-lipid LNP
 - Many different lipids:
 - Phospholipids (DOPE, DPyPE, DSPC, DPPC, DOPC, PLPE)
 - PEGylated lipids (DMG-PEG2k, DSG-PEG2k, ALC-0159, DSPE-PEG2k, DSPE-PEG5k)
 - Sterol lipids (Cholesterol, β -sitosterol)
 - Ionizable lipids (DODMA, Dlin-MC3-DMA, SM-102, ALC-0315)
 - Different nucleic acids:
 - mRNA, saRNA, gRNA, pDNA, etc.




IM21.7c is suitable with Different Phospholipids for LNP formulation

IM21.7c (mM)	DODMA (mM)	Helper Lipid (mM)	Cholesterol (mM)	DSG-PEG _{2k} (mM)	Size	PDI	Zeta	EE%
IM21.7c 4	DODMA 3	DPyPE	1	0.15	54±1	0.085	+16	100
		PLPE	1		52±9	0.197	+14	100
		DOPE	1		60±5	0.155	+13	100
		DOPC	1		50±10	0.163	+17	100
		DPPC	1		59±8	0.134	+17	100
		DSPC	1		43±3	0.162	+15	100



5-lipid LNP formulation


Caco-2 cells – Luciferase expression (500 ng)



IM21.7c is suitable with Different PEGylated Lipids for LNP formulation

Cationic lipid (mM)	Ionizable lipid (mM)	Helper Lipid (mM)	Cholesterol (mM)	PEG-lipid (mM)		Size	PDI	Zeta	EE%
IM21.7c 4	DODMA 3	DPyPE 1	1.85	DSG-PEG _{2k}	0.15	75±5	0.050	+11	100
				DSPE-PEG _{5k}	0.15	61±4	0.082	+2	100
				DSPE-PEG _{2k}	0.15	85±3	0.056	+24	100
				DMG-PEG _{2k}	0.15	40±3	0.198	+18	100
				ALC-0159	0.15	36±6	0.216	+18	100


Caco-2 cells – Luciferase expression (500 ng)



5-lipid LNP formulation



IM21.7c is suitable with different Sterols for LNP formulation

Cationic Lipid (mM)	Ionizable Lipid (mM)	Helper Lipid (mM)	Sterol (mM)	DSG-PEG _{2k} (mM)	Size	PDI	Zeta	EE%
IM21.7c 4	DODMA 3	DPyPE 1	Cholesterol	1.85	55±3	0.152	+14	100
			β-sitosterol	1.85	77±9	0.163	+14	100
			Stigmasterol	1.85	127±11	0.177	+7	100




Cholesterol

5-lipid LNP formulation



β-sitosterol



Stigmasterol

Screening mRNA-LNP formulations

Table 1: Summary of CQAs for mRNA product development and manufacture (continued)

Attribute quality	CQA*	Justification of importance	Stage tested	Stage of manufacture (DS/DP)	Is CQA unique for RNA products?	Regulatory precedence?	Release and stability assay methods
Product characterization and quality	Encapsulation efficiency	Encapsulation protects the mRNA from degradation and stabilizes it, it is needed for delivery and efficacy, and is an important measure of quality	Stability/Release	DP	No	Yes	RP-LC-Fluorescence detection
	Surface charge	Surface charge affects the distribution of nanoparticles and affects performance and process consistency	Characterization	DP	No	Yes	ELS iCIEF TNS binding assay
	LNP size	Affects distribution and uptake of nanoparticles and impacts safety and efficacy	Stability/Release	DP	No	Yes	DLS (SR-DLS or RT-MALS in/online option) NTA
	LNP polydispersity	Polydispersity is a measure of the homogeneity and impacts safety and efficacy	Stability/Release	DP	No	Yes	DLS (SR-DLS or RT-MALS in/online option)

<https://www.biophorum.com/news/an-industry-standard-for-mrna-lnp-analytics/>

Screening mRNA-LNP formulations

Table 1: Summary of CQAs for mRNA product development and manufacture (continued)

Attribute quality	CQA*	Justification of importance	Stage tested	Stage of manufacture (DS/DP)	Is CQA unique for RNA products?	Regulatory precedence?	Release and stability assay methods
Product characterization and quality	Encapsulation efficiency	Encapsulation protects the mRNA from degradation and stabilizes it, it is needed for delivery and efficacy, and is an important measure of quality	Stability/Release	DP	No	Yes	RP-LC-Fluorescence detection
	Surface charge	Surface charge affects the distribution of nanoparticles and affects performance and process consistency	Characterization	DP	No	Yes	ELS iCIEF TNS binding assay
	LNP size	Affects distribution and uptake of nanoparticles and impacts safety and efficacy	Stability/Release	DP	No	Yes	DLS (SR-DLS or RT-MALS in/online option) NTA
	LNP polydispersity	Polydispersity is a measure of the homogeneity and impacts safety and efficacy	Stability/Release	DP	No	Yes	DLS (SR-DLS or RT-MALS in/online option)

EE%: > 95%

Zeta potential: 5-20 mV

Size: 25-100 nm

PDI: < 0.2

<https://www.biophorum.com/news/an-industry-standard-for-mrna-lnp-analytics/>

LipidBrick® library: mRNA-LNP characterization

Cationic / Ionizable lipid (mM)		Helper Lipid (mM)		Cholesterol (mM)	DSG-PEG (mM)
IMXX.Xc/DODMA		4/3	DPyPE	1	1.85
					0.15

	IM21.7c	IM16c	IM13c	IM15c	IM22c	IM12c	IM25c	IM3c
Size (nm)	45 ± 2	56 ± 3	74 ± 4	66 ± 6	40 ± 2	29 ± 1	86 ± 2	41 ± 2
Zeta (mV)	+10	+13	+11	+9	+14	+12	+10	+16
PDI	0.158	0.132	0.090	0.113	0.140	0.120	0.053	0.104
EE (%)	99	99	99	99	99	99	99	99

Results confirmed in 4-lipid LNP formulation.

Screening mRNA-LNP formulations

Table 1: Summary of CQAs for mRNA product development and manufacture (continued)

Attribute quality	CQA*	Justification of importance	Stage tested	Stage of manufacture (DS/DP)	Is CQA unique for RNA products?	Regulatory precedence?	Release and stability assay methods
Product characterization and quality	Encapsulation efficiency	Encapsulation protects the mRNA from degradation and stabilizes it, it is needed for delivery and efficacy, and is an important measure of quality	Stability/Release	DP	No	Yes	RP-LC-Fluorescence detection
	Surface charge	Surface charge affects the distribution of nanoparticles and affects performance and process consistency	Characterization	DP	No	Yes	ELS iCIEF TNS binding assay
	LNP size	Affects distribution and uptake of nanoparticles and impacts safety and efficacy	Stability/Release	DP	No	Yes	DLS (SR-DLS or RT-MALS in/online option) NTA
	LNP polydispersity	Polydispersity is a measure of the homogeneity and impacts safety and efficacy	Stability/Release	DP	No	Yes	DLS (SR-DLS or RT-MALS in/online option)
	Lipid content	Lipid content impacts function of the DP	Stability/Release	DP	No	Yes	LC-CAD LC-ELSD

Quantification by BIA experts

<https://www.biophorum.com/news/an-industry-standard-for-mrna-lnp-analytics/>

Lipid quantification (preliminary results)

	Theoretical mass%	Measured mass%	Starting mass	Mass yield
IM21.7c	63.0%	66.5%	1.802 mg	1.515 mg
Cholesterol	21.1%	18.2%	0.603 mg	0.418 mg
DOPE	10.5%	11.0%	0.301 mg	0.257 mg
PEG-DMG 2K	5.3%	4.3%	0.152 mg	0.099 mg

Data generated by Sartorius BIA Separations.


Analysis was done on a PATfix™ system with an ELSD detector, using a monolithic CIMac™ C4 HLD-0.1 Analytical Column (2 µm).

Mass yield is 84%, not a lot of lipids are lost.


Ratio are slightly different in favour of the cationic lipid and the DOPE.

LipidBrick® enables to formulate mRNA-LNPs


4-lipids LNP IM21.7c



5-lipids LNP IM21.7c



SM-102 LNP (Moderna)




CryoTEM images

Get Beyond the Liver with LipidBrick®


LipidBrick® Library Offers Consistent Efficiency & Flexibility

- The sole modification of the cationic lipid leads to a different efficiency
- Repeatable across different cell lines




LipidBrick® Library Offers Consistent Efficiency & Flexibility



- Depending on the cell type, transfection efficiency differs simply by changing the cationic lipid of the same LNP formulation
- LipidBrick® library gives more options




Adapt LNP Tropism to Your Specific Application with LipidBrick® Library




- Deliver RNA to different organs while **bypassing liver**
- Enhance the efficacy of LNPs by **fine-tuning formulations** around a specific LipidBrick® (here IM21.7c)
- LipidBrick® Library provides new delivery options for RNA-based therapeutics & vaccines



Get Beyond the Liver with LipidBrick® (in 4- or 5-lipid LNP conformation)



IM21.7c-LNPs mainly targets endothelial cells and alveolar macrophages



Unpaired t test *p<0.5 and **p<0.01


Reduce Liver Off-Target Effects with LipidBrick®

Intramuscular injection



- Precise local targeting
- Reduced Side Effects
- Immune Resp. Management
- Interesting for oncology, vaccination, local gene editing or protein replacement, tissue regeneration

Luciferase expression after IM injection (24 h - 5 µg/mouse)




Transfect the Lungs Efficiently using LipidBrick® via Intranasal Instillation

Intranasal instillation





- Interesting for treatments for cystic fibrosis



LipidBrick[®]-LNPs are safe &
weakly immunogenic


IM21.7c-LNPs triggers low activations of primary human moDCs *in vitro*



LipidBrick®-LNPs triggers the same low cytokine profile as ALC-0315 LNPs




IV injection
CleanCap® Fluc mRNA (5moU); 10µg



Cationic/Ionizable lipid (mM)		Helper lipid (mM)		Cholesterol (mM)	PEG lipid (mM)	
LNP LipidBrick® 5-lipid (IMXXc/DODMA)	4/3	DPyPE	1	1.85	DSG-PEG _{2k}	0.15

Average (N=6 mice)			PBS			LPS			LNP BNT			LNP with LipidBrick® IM21.7c			LNP with LipidBrick® IM12			Legend in pg/mL
Function	Method	Cytokine	2 h	6 h	24 h	2 h	6 h	24 h	2 h	6 h	24 h	2 h	6 h	24 h	2 h	6 h	24 h	
Pro-inflammation Innate & Adaptive immunity Anti-viral and anti-tumoral activities Cytokine release	MACSPlex	IFN-γ																> 10 ⁴
Pro-inflammation Cytokine release Cell proliferation Cell Proliferation & Apoptosis Anti-viral and anti-tumoral activities	MACSPlex	TNF-α																10 ³ - 10 ⁴
Pro-inflammation Macrophage activation Increase neutrophil & monocyte function	MACSPlex	GM-CSF																201-1000
Pro-inflammation	ELISA	IL-1α																51-200
Pro-inflammation Cell Proliferation & Apoptosis	ELISA	IL-1β																0-50
Regulate T-cell proliferation and differentiation	MACSPlex	IL-2																
Anti-inflammation T-cell & B-cell proliferation & differentiation	MACSPlex	IL-4																
Pro-inflammation B-cell differentiation Eosinophil activation	MACSPlex	IL-5																
Pro-inflammation Differentiation Cytokine release	ELISA	IL-6																
Anti-inflammation Inhibition of pro-inflammatory cytokines Immunoregulation	MACSPlex	IL-10																
Pro-inflammation Cell differentiation Activates NK cells Regulates T-cell & B-cell cooperation	MACSPlex	IL-12																
Pro-inflammation	MACSPlex	IL-17A																
Pro-inflammation	MACSPlex	IL-23																


LipidBrick®-LNPs triggers the same low cytokine profile as ALC-0315 LNPs




Average (N=6)			LNP BioNTech ALC-0315			LNP LipidBrick® IM21.7c 5-lipid			LNP LipidBrick® IM21.7c 4-lipid			LNP LipidBrick® IM12c 5-lipid			Legend in pg/ml
Function	Method	Cytokine	2h	6h	24h	2h	6h	24h	2h	6h	24h	2h	6h	24h	
Pro-inflammation Innate & Adaptative immunity Anti-viral and anti-tumoral activities Cytokine release	MACSPlex	IFN-γ													> 10 ⁴
Pro-inflammation Cytokine release Cell proliferation Cell Proliferation & Apoptosis Anti-viral and anti-tumoral activities															
Pro-inflammation Macrophage activation Increase neutrophil & monocyte function															
Pro-inflammation															
Pro-inflammation Cell Proliferation & Apoptosis	ELISA	TNF-α													10 ³ - 10 ⁴
Regulate T-cell proliferation and differentiation															
Anti-inflammation T-cell & B-cell proliferation & differentiation	MACSPlex	GM-CSF													201-1000
Pro-inflammation B-cell differentiation Eosinophil activation															
Pro-inflammation Differentiation Cytokine release	ELISA	IL-1α													51-200
Anti-inflammation Inhibition of pro-inflammatory cytokines Immunoregulation															
Pro-inflammation Cell differentiation Activates NK cells Regulates T-cell & B-cell cooperation	MACSPlex	IL-1β													0-50
Pro-inflammation															
Pro-inflammation T-cell & B-cell proliferation & differentiation	MACSPlex	IL-2													
Pro-inflammation B-cell differentiation Eosinophil activation															
Pro-inflammation Differentiation Cytokine release	ELISA	IL-4													
Anti-inflammation Inhibition of pro-inflammatory cytokines Immunoregulation															
Pro-inflammation Cell differentiation Activates NK cells Regulates T-cell & B-cell cooperation	MACSPlex	IL-5													
Pro-inflammation															
Pro-inflammation Differentiation Cytokine release	ELISA	IL-6													
Anti-inflammation Inhibition of pro-inflammatory cytokines Immunoregulation															
Pro-inflammation Cell differentiation Activates NK cells Regulates T-cell & B-cell cooperation	MACSPlex	IL-10													
Pro-inflammation															
Pro-inflammation	MACSPlex	IL-12													
Pro-inflammation															
Pro-inflammation	MACSPlex	IL-17A													
Pro-inflammation															
Pro-inflammation	MACSPlex	IL-23													
Pro-inflammation															

LipidBrick® for Vaccination & Immunisation

IM21.7c-based LNPs lead to strong humoral immune response




IM12c-based LNPs lead to strong humoral immune response




LipidBrick® for *Ex Vivo* Cell Therapy

Reach High Ex Vivo Delivery Efficacy & Viability with IM21.7c & IM25c-LNPs on Human Primary T cells



Results obtained after 48 hours; LipidBrick-LNPs do not require ApoE supplementation
Dedicated 4-lipid LNP formulation + Optimized Ex vivo T cells protocol available

IM21.7c-LNPs lead to Efficient DNA Transfection of Primary T cells




Results obtained after 48 hours; LipidBrick-LNPs do not require ApoE supplementation
Dedicated 4-lipid LNP formulation + Optimized Ex vivo T cells protocol

Cationic/Ionizable lipid (mM)		Helper Lipid (mM)		Cholesterol (mM)	PEG lipid (mM)	
LNP IM21.7c 4-lipid		2	DOPE	3	4.85	DMG-PEG _{2k} 0.15

IM21.7c-LNPs are Stable


IM21.7c-LNPs maintain Size, PDI & zeta potential stable for ≥3 months at 4°C & -20°C



4-lipid LNP formulation
Buffer containing sucrose

Cationic/Ionizable lipid (mM)		Helper Lipid (mM)		Cholesterol (mM)	PEG lipid (mM)	
LNP Onpatto (Dlin-MC3)	5	DSPC	1	3.85	DMG-PEG _{2k}	0.15
LNP LipidBrick® IM21.7c	5	DOPE	1	3.85	DMG-PEG _{2k}	0.15

IM21.7c-LNPs Maintain a stable EE% & Transfection Efficiency for ≥3 months at -20°C



4-lipid LNP formulation
Buffer containing sucrose

Cationic/Ionizable lipid (mM)		Helper Lipid (mM)		Cholesterol (mM)	PEG lipid (mM)	
LNP Onpatto (Dlin-MC3)		5	DSPC	1	3.85	DMG-PEG _{2k} 0.15
LNP LipidBrick® IM21.7c		5	DOPE	1	3.85	DMG-PEG _{2k} 0.15

Summary

LipidBrick® advantages

- **Efficient:** Modulate LNP properties to adapt the biodistribution depending on the therapeutic purpose
 - LipidBrick® shows a great potential by adapting and enlarging the biodistribution of LNPs
- **Secure:** Use a unique lipid structure protected by an independent patent owned by Polyplus®
 - Independent IP covering the use of LipidBrick® library with nucleic acid in LNP formulation
- **Support:** *in vitro* and *in vivo* proof of concept have been successfully performed
 - In addition to our extensive PoC data, our Scientific Support team is here to help you optimize and accelerate the development of your LNPs.
- **Flexible:** Versatile Library allowing to select the Best Suited LipidBrick® for each application
 - Consistent Efficiency & Flexibility to Fine-Tune your LNP formulations on your specific therapeutic purposes

LipidBrick® Library kit (8 x 50mg) is now available

- **Flexibility:** Select from 8 unique, IP-protected LipidBrick® structures for complete application coverage
- **Consistency:** All shown to be efficient in LNP formulations for *in vitro* & *in vivo* applications
- **Optimization:** Find your ideal LipidBrick® fast for each of your therapeutic application
- **Perfect trial size:** Accessible, convenient, and ideal for 1st screening rounds



Item	Description	Type	Quality Grade	Pack size	Order number
LipidBrick® Library Kit	50 mg of each of 8 different LipidBrick® lipids	Cationic lipids (powder and oil)	Research Grade*	8 lipids x 50 mg	101000241

*GMP grade development of specific lipid(s) upon request

Accessibility Leads to Accelerated Therapeutics

- Earlier access to lipids that can be clinically translated allows for accelerated development of therapeutics
- Sartorius is committed to flexible purchase & use agreements
 - Use under standard purchase and supply agreements
 - Flexible licensing agreements
 - Flexible Development Options for library lipids
- Removing the financial barrier to LNP formulation development

Product	Part number	Pack size
LipidBrick® IM21.7c	101000232	50 mg - NEW
	101000172	250 mg
	101000173	1 g
LipidBrick® Library Kit	101000241	8 lipids x 50 mg - NEW



Removing the Barrier to Entry- Formulation Support Built In

- **Expert Scientific & Regulatory Support teams**

- Expertise in LNP Formulation, High-Throughput Screening (HTS), Design of Experiments (DoE) & Transfection
- Access to a Dedicated Team of Field Application Specialists for Specific On-Site Support
- Scientific, Quality and Regulatory assistance (>15 years of GMP Expertise in Excipients & Raw Materials)

- **Optimized LNP Formulations Guidance**

- Recommended 'starting points' LNP formulations for enhanced performance in various applications (*in vivo*, *ex vivo*, etc.)
- Backed by a wealth of internal Proof of Concept data


- **Protocols & Methodologies**

- Detailed protocols for LNP Formulation and Characterization (DLS)
- Step-by-step Transfection Guidelines for *in vitro*, *in vivo*, and *ex vivo* applications

- **In-Lab Formulation Tools**

- Access to our LNP formulation spreadsheet to Ease the Generation of new LNP formulations (and associated calculations)

Any questions?



Quentin Huaulmé, PhD
Chemistry Research Scientist
Polyplus, now part of Sartorius
Quentin.huaulme@sartorius.com

Claire Guéguen, PhD
R&D Life Sciences Manager
Polyplus, now part of Sartorius
Claire.gueguen@sartorius.com

SARTORIUS