


Elucidation of Lipid Nanoparticle Fusion Mechanisms Using Artificial Membranes

Shu Okumura, Kozue Furuya, Mitsuko Sugano

Toshiba Corporation, Corporate Laboratory, Materials & Frontier Research Center, Frontier Research Department

Lipid Nanoparticle (LNP) Problems and “Designable LNP” Concept

Cell Target LNP



Problem

Targeted delivery to specific cell types is essential for advancing LNP clinical applications.

Conventional Solution: Active Targeting

Conjugate specific antibodies or ligands to LNPs




Suppress immunogenic responses and enhance application versatility.

Our Solution


We engineered “Designable LNP” as an LNP design platform to achieve cell-tropism by **optimizing lipid composition**.

"Designable LNP" Platform

Optimized design method for LNP composition




Design of Cell-tropic LNP



Why LNPs have cell tropism?

Cells exhibit distinct lipid compositions



- AI-driven lipid composition optimization **cannot** identify the underlying causes of cell specificity.

• **Research Question:** Investigate how the cell-tropic delivery efficiency of LNPs correlates with the lipid composition of each cell type.

Lipid-lipid interaction

distinct lipid profiles of each cell type can be separated.


Objective

- Elucidate the cell-tropic delivery mechanism of "Designable LNP".
- Verify the cell-tropic delivery mechanism based on **lipid-lipid interactions**, excluding non-lipid membrane components.
- Construct a cell-tropic delivery analysis using **cell-mimic nanoparticles**.

Materials


Lipid Nanoparticles

LNPs




Ionized lipids

FFT-10



FFT-20




Chol = Cholesterol PEG = DMG-PEG-2000

	LNP1	LNP2	LNP3	LNP4	LNP5	LNP6
FFT10	26%	15%	41%	28%	22%	0%
FFT20	0%	0%	0%	0%	0%	27%
DOPE	0%	0%	0%	0%	0%	5%
DOTAP	18%	21%	15%	10%	31%	10%
Chol	53%	60%	42%	58%	44%	55%
PEG	4%	4%	3%	4%	3%	4%
z-ave(nm)	89.8	70.44	111.7	78.78	84.89	71.02
ζ(mV)	19.0	13.3	21.6	4.71	19.9	11.1

Cell mimic nanoparticles

SUVs




	A	B	C	D	E	F	G	H	I
DOPC	63%	53%	43%	58%	48%	38%	53%	43%	33%
DOPE	0%	0%	0%	5%	5%	5%	10%	10%	10%
DOPS	5%	5%	5%	5%	5%	10%	10%	10%	10%
Chol	30%	40%	50%	30%	40%	50%	30%	40%	50%
PEG	1%	1%	1%	1%	1%	1%	1%	1%	1%

	J	K	L	M	N	O	P	Q	R
DOPC	58%	48%	38%	53%	43%	33%	48%	38%	28%
DOPE	0%	0%	0%	5%	5%	5%	10%	10%	10%
DOPS	10%	10%	20%	20%	20%	20%	20%	30%	30%
Chol	30%	40%	50%	30%	40%	50%	30%	40%	50%
PEG	1%	1%	1%	1%	1%	1%	1%	1%	1%


Methods

Mixing NPs

LNP pattern



SUV pattern



Fusion

Dynamic Light Scattering

z(nm) average Size measurement
By Dynamic Light Scattering

Particle size UP

Measured by Zetasizer
(Malvern. Inc.)

Statistical Analysis

$$R = \frac{D_{\text{mix}}}{\sqrt{D_{\text{LNP}}^2 + D_{\text{SUV}}^2}}$$

- D_{LNP} : Diameter of the liposome before mixing
- D_{SUV} : Diameter of the small unilamellar vesicle (SUV) before mixing
- D_{mix} : Diameter of the particle after mixing
- R : Particle size change rate


Model1: SUV composition +
LNP composition model

Model2: SUV composition +
LNP composition + LNP-SUV
interaction model

Characteristics of mGFP LNP

RNA-LNP uptake in cells

Scale Bar = 100 μ m




LC-MSMS Cell lipid compositions analysis

LC-MSMS Cell lipid compositions analysis

Lipids(mol %)	Jurkat	MDA-MB-231	Huh-7
PC	48.8	50.1	48.9
PS	1.1	5.1	2.3
PE	2.9	13.4	6.1
PI	2.7	2.3	7.9
SM	17.2	6.8	8.2
Cholesterol	27.2	22.3	26.6
Total	100.0	100.0	100.0

RNA-LNP uptake in cells



LNP1 ⇒ Standard

LNP2,4 ⇒ Lower uptake


LNP3 ⇒ High uptake in All types

LNP5 ⇒ High uptake in **Huh7** ⇒ High anionic/cholesterol?


LNP6 ⇒ High uptake in **Jurkat** ⇒ Low anionic?

Statistical Analysis of Particle Size Change Rate Using Linear Regression

Model 1: SUV Composition + LNP Composition



Model 2: SUV Composition + LNP Composition + LNP-SUV Interaction



Introducing the **LNP-SUV interaction** led to improved performance.

Conclusion and Future Prospects

Conclusion

- We confirmed a cell-tropic delivery mechanism based on **lipid-lipid interactions**, excluding the influence of non-lipid molecules.
- We constructed a cell-tropic delivery analysis model using cell-mimic nanoparticles.



Future Prospects

- Employing various nanoparticles, we will comprehensively explore factors involved in cell-tropism.



APPENDIX

Statistical Analysis of Particle Size Change Rate Using Linear Regression



(1) SUV composition + LNP composition model





(2) SUV composition + LNP composition + intracellular interaction model



(3) SUV composition + LNP composition + LNP-SUV interaction model




(4) SUV composition + LNP composition + LNP-SUV interaction + intracellular interaction model



Property of LNPs

	z (nm)	Pdl	ζ potential(mV)
LNP1	89.8	0.145	19
LNP2	70.44	0.229	13.3
LNP3	111.7	0.106	21.6
LNP4	78.78	0.148	4.71
LNP5	84.89	0.168	19.9
LNP6	71.02	0.133	11.1

Characteristics of LNP



Scale Bar = 100 μ m