


ULTRA-LONG-ACTING BIODEGRADABLE POLYMERIC SOLID IMPLANTS FOR HIV TREATMENT MAINTENANCE

Annu Anna Thomas, MS, Ph.D.
Benhabbour Research Laboratory
Lampe Joint Department of Biomedical Engineering
The University of North Carolina at Chapel Hill


Hurdles In HIV Treatment Progress

Problem

- 39 million people are living with HIV.
- 53% of all people living with HIV are women and girls.
- Heavy pill burden leading to poor patient adherence.





39.9 million
people living with HIV
[36.1-44.6 million]



Figures created using Bio render, adapted from source **WHO, 2024** & **NIAID**

Bridging Gaps in HIV Treatment Maintenance



Cabenuva (Cabotegravir 600 mg/3mL & Rilpivirine 900 mg/3 mL) is not removable, up to 6 injections a year, painful IM injections, missed doses can lead to drug resistance.


Figures created using Bio render, Source: <https://www.cabenuva.com/>

Bridging Gaps in HIV Treatment Maintenance

Solution

Polymeric Solid Implants (PSI)

- Technology derived from *Injectable in-situ forming implants (ISFI)*.



Advantages:


1. Ultra-Long-Acting ($\geq 1\text{year}$).
2. Solvent-free (low toxicity), biodegradable and removal if treatment termination is required.
3. High drug loading($\sim 90\text{ wt\%}$), accommodate multiple drugs, variability in size and shape.

Figures created using Bio render, Source: Maturavongsadit, Panita, et al. "Biodegradable polymeric solid implants for ultra-long-acting delivery of single or multiple antiretroviral drugs." *International journal of pharmaceutics* 605 (2021): 120844. Benhabbour, S. Rahima, et al. "Ultra-long-acting tunable biodegradable and removable controlled release implants for drug delivery." *Nature communications* 10.1 (2019): 4324.

Antiretroviral loaded PSIs for HIV Treatment Maintenance


Overcome current limitations of Long-Acting therapies

- Cabotegravir and Rilpivirine loaded PSIs.
- PSIs API loading comparable to Cabenuva.
- Achieve *in vitro* release for longer duration (> 6 months) to minimize clinic visits and enhance adherence.
- PSIs capable of delivering CAB and RPV *in vivo* above their targeted plasma drug concentration (> 4 × PA-IC₉₀).



Cabotegravir and Rilpivirine PSIs for HIV Treatment Maintenance

Targeted Drug loading




Physical Dimensions and Weights of CAB and RPV PSIs

PSI	Analytical Drug loading (mg/g)	Weight (mg)	Length (mm)	Height (mm)	Width (mm)	Surface Area (mm ²)	Specific Surface Area (mm ⁻¹)
CAB PSIs	886.76± 0.79	102.70	10	3.01	2.80	132.80	1.58
		106.81	10	3.25	2.74	137.61	1.54
		99.60	10	2.99	2.86	134.10	1.56
		107.91	10	3.70	2.81	151.01	1.45
RPV PSIs	900.01± 17.22	96.20	10	3.68	2.80	150.20	1.46
		102.90	10	3.69	2.81	150.74	1.45

Cabotegravir and Rilpivirine PSIs for HIV Treatment Maintenance

Targeted Drug loading and in vitro release study



Drug	Drug loading/PSI (wt%)	Total drug per PSI (mg)	Burst release at 24 h (%)	Burst release at 24 h (µg)	Cumulative µg/day release (1-7 d)	Cumulative µg/day release (7-60 d)	Cumulative µg/day release (60-365 d)	Estimated time to 100% release (days)
CAB	90	91.36 ± 3.2	0.17 ± 0.06	157.38 ± 63.42	87.31	64.88	61.22	~ 1918 days
RPV	90	92.09 ± 5.28	0.03 ± 0.01	27.48 ± 6.36	20.75	25.67	35.71	~ 3134 days

Cabotegravir and Rilpivirine PSIs for HIV Treatment Maintenance

In vivo PK Study




Study Plan

Treatments (CAB/RPV PSI SC Implants):

- G.1- CAB PSI: 4500 mg/kg
- G.2- CAB PSI: 2250 mg/kg
- G.3- RPV PSI: 4500 mg/kg
- G.4- RPV PSI: 2250 mg/kg
- G.5- CAB + RPV PSI: 4500 mg/kg each


Plasma collection to quantify drug concentration (PK)

Collection of residual implant to quantify remaining drugs




Cabotegravir and Rilpivirine PSIs for HIV Treatment Maintenance

In vivo PK Study 180 d



**PSIs retrieved at d180
post administration**




Group	API	Dose (mg/kg)	Administered PSI average weight (mg)	~ % residual drug	~ % PLGA Degradation
G1	CAB	4500	100.21 ± 2.57	85.66 ± 0.19	87.25
G2	CAB	2250	100.70 ± 1.50	78.32 ± 5.65	90.81
G3	RPV	4500	99.78 ± 1.21	83.43 ± 10.79	88.94
G4	RPV	2250	103.30 ± 1.27	82.09 ± 2.49	88.75
G5	CAB	4500	104.11 ± 2.43	94.41 ± 3.83	91.37
G5	RPV	4500	103.70 ± 1.52	90.86 ± 2.26	90.69

*N=2-3 for each group

Conclusion

- Developed Polymeric Solid Implants that can accommodate antiretrovirals (CAB ~900 mg/g, and RPV ~900 mg/g).
- Sustained zero order release kinetics *in vitro* for ~365 days.
- CAB/RPV PSIs elicited sustained CAB and RPV release with plasma drug concentrations $>4 \times \text{PA-IC}_{90}$ (664 ng/ml – CAB; 48 ng/ml – RPV) for 180 days.



Acknowledgements

Benhabbour Lab

Dr. S. Rahima Benhabbour
Dr. Nanditha C K
Dr. Jasmine King
Thy Le
Caroline Dau
Kayla Marie Maue

UNC Collaborators

Dr. Mackenzie Cottrell
Dr. Craig Sykes
Amanda Schauer



1R01AI176949-01A1

Lampe Joint Department of Biomedical Engineering



PHILADELPHIA, PA
JULY 13-18, 2025

**NEXT-GENERATION
DELIVERY INNOVATIONS**