

Tech Session 2: Delivery Technologies in Consumer and Diversified Industrial Products(C&DP 1)

FATIMA ABID

CRS 2022 Annual Meeting & Expo

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada

Advanced Delivery Science



University of
South Australia

Pharmaceutical Innovation and Development Group

Your Partner in Product Development Research





University of
South Australia

Antibiotic discovery for Resistant Acne Infections: Novel topical formulation development

F. Abid¹, A. Parikh¹, Y. Song¹, S. H Youssef¹, S. W Page², DJ. Trott³ and S. Garg¹

¹ University of South Australia, Adelaide, South Australia, 2042, Australia

² Advanced Veterinary Therapeutics, Newtown, New South Wales, 3125, Australia

³ The University of Adelaide, Roseworthy, South Australia, 5371, Australia

OUTLINE

- ❖ Introduction of model drug candidate ADS1
- ❖ Formulation development
- ❖ Physicochemical characterisation
- ❖ In vitro evaluation of encapsulated ADS1
- ❖ Conclusion and future directions
- ❖ Acknowledgement



❖ Exhibit Antimicrobial properties

- Anti-viral
- Anti-fungal
- Antibacterial

❖ Anticancer properties



REPOSITIONING

Cutibacterium acnes (C. acnes)


$\mu\text{g/mL}$	SDA321
MIC_{50}	0.031
MIC_{90}	> 0.062
MIC Range	$\leq 0.008 - 0.062$

Table 1. MIC_{50} , MIC_{90} and MIC range against *C. acnes* isolates

Staph aureus (S. aureus) and Staph epidermidis (S. epidermidis)

$\mu\text{g/mL}$	SDA321
<i>S. aureus</i> ATCC 29213	0.25
<i>S. epidermidis</i> ATCC 14990	0.25

Table 2. MIC range against *S. aureus* and *S. epidermidis* isolates



ANALYTICAL METHOD DEVELOPMENT

ADS1 Lacks chromophore




Fig 1: ADS1 Chromatogram (700 µg/mL)

Column	Eclipse Plus C18 (Agilent), 5µm, 4.6 x 150mm
Mobile Phase	Methanol: Water: Glacial Acetic Acid (92:8:1)
Flow rate	1 mL/min
Injection volume	50 µl
Temperature	30 C

Table 3: Method development parameters

Parameters	ADS1
Linearity range	200 – 900 µg/mL
Correlation coefficient (r)	1.0
Accuracy (Mean \pm SD)	99.98 % \pm 1.02
LOD	7.1 µg/ml
LOQ	21.50 µg/ml
Precision	< 2% (Pass)
Repeatability/intermediate precision	RSD < 2% (Pass)
Robustness (Inj vol, temp)	RSD < 2% (Pass)
Systems suitability	Compliant

Table 4: Method validation parameters

GREENNESS EVALUATION




Fig 2: AGREE and GAPI pictogram of reported HPLC method




Fig 3: AGREE and GAPI pictogram of reported LCMS method





Fig 4: AGREE and GAPI pictogram of developed HPLC method



CONSIDERATIONS FOR FORMULATION DEVELOPMENT


Poor aqueous solubility ($<40 \mu\text{g/mL}$)

Nano-sized carriers for efficient hair follicular delivery (SNMSD)



FORMULATIO

N DEVELOPMEN T



- ❖ Polymer selection (Soluplus, HPMC (ASLG, ASMG)
- ❖ Drug: Polymer ratio (1:5)

F1	ADS1-Soluplus
F2	ADS1-HPMC ASLG
F3	ADS1-HPMC ASMG



SOLUBILITY ENHANCEMENT




Fig 5. Solubility enhancement study



PHYSICOCHEMICAL CHARACTERISATION (1)




Fig 6. TEM image




Fig 7. Particle size using DLS (71 nm)






Fig 8. SEM study (A: ADS1; B: Soluplus; C: F1)

PHYSICOCHEMICAL CHARACTERISATION (2)




Fig 9. DSC study

Pharmaceutical Innovation and Development Group

Your Partner in Product Development Research



University of
South Australia



PHYSICOCHEMICAL CHARACTERISATION (3)




Fig 10. FTIR spectra (A: F1; B: Physical mix (ADS1-Soluplus); C: Soluplus; D: ADS1)



In vitro evaluation: Antibacterial activity and cell viability




Table 5. MIC₅₀, MIC₉₀ and MIC against *Cutibacterium acnes* (*C. acnes*)

µg/mL	ADS1	F1	F2	F3
MIC ₅₀	0.031	0.031	0.062	0.062
MIC ₉₀	> 0.062	> 0.062	> 0.062	> 0.062
MIC Range	≤ 0.008 – 0.062	≤ 0.008 – 0.062	≤ 0.008 – 0.062	≤ 0.008 – 0.062

Table 6. MIC range against *S. aureus* and *S. epidermidis*

µg/mL	ADS1	F1	F2	F3
<i>S. aureus</i> ATCC 29213	0.25	0.5	0.5	0.5
<i>S. epidermidis</i> ATCC 14990	0.25	0.5	0.5	0.5





Fig 11. Percentage viability (A:HFF-1 cell line, B: HaCaT cell line)



CONCLUSION


- ❖ ADS1 is an effective antibacterial agent
- ❖ Self-nanomicellising solid dispersion presents a promising drug delivery technology for poorly water-soluble drugs
- ❖ Findings from in vitro evaluations in both bacterial and cellular systems warrant further studies with potential transfer into suitable in vivo models.
- ❖ Ongoing work: optimisation and further characterisation of ADS1-SNMSD incorporating transparent gel



ACKNOWLEDGEMENT



University of
South Australia



Prof Sanjay Garg
Dr Yunmei Song
Dr Ankit Parikh

Future Industries Institute

Prof Darren Trott

Adelaide Microscopy

Dr Stephen Page



Pharmaceutical Innovation and Development Group



CRS 2022 Annual Meeting & Expo

Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada



THANK YOU!



Fatima.abid@mymail.unisa.edu.au



@PIDG_UniSA



@FatimaAbid91

<https://www.linkedin.com/in/fatima-bd/>



CRS 2022 Annual Meeting & Expo

Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada