

Enabling the Rational Development of Long-Acting Contraceptive Levonorgestrel Intrauterine Systems

SURAJ FANSE

UNIVERSITY OF CONNECTICUT

CONTROLLED RELEASE SOCIETY


CRS 2023 ANNUAL MEETING & EXPOSITION

JULY 24-28, 2023

Paris Hotel » **Las Vegas, NV, USA**

THE FUTURE OF DELIVERY SCIENCE

Levonorgestrel intrauterine systems



Levonorgestrel intrauterine systems (LNG-IUSs)
Long-Acting Reversible Contraceptives (LARCs)

Hormonal intrauterine systems for contraception

Advantages

- ✓ >99% effective
- ✓ Reversible
- ✓ Non-surgical
- ✓ Ultra long-acting (3-8 years)

Commercial LNG-IUS products



Product name	Strength	Effective duration	<i>In vivo</i> drug release rate at the end of first year of use	<i>In vivo</i> drug release rate at the end of approved year of use	Applicant holder	Year of approval
Mirena®	52 mg	7 years*	18.0 µg/ day	8.0 µg/ day (7 years)	Bayer	2000
Liletta®	52 mg	6 years	17.0 µg/ day	8.6 µg/ day (6 years)	Allergan/ Medicines360	2015
Kyleena®	19.5 mg	5 years	9.8 µg / day	7.4 µg/ day (5 years)	Bayer	2016
Skyla®	13.5 mg	3 years	6.0 µg/ day	5 µg/ day (3 years)	Bayer	2013

*The US FDA has extended the duration of use for Mirena® up to 8 years in 2022.

Generic drug product development of LNG-IUSs is challenging!!

Challenges

Why is generic product development of LNG-IUSs challenging?



Objectives of current research

Impact of critical material attributes on product performance

Influence of processing parameters on product performance

Elucidating drug release mechanisms

Establishing *in vitro* bioequivalence testing methods




CRS 2023 ANNUAL MEETING & EXPOSITION


JULY 24-28, 2023 **Paris Hotel » Las Vegas, NV, USA**

THE FUTURE OF DELIVERY SCIENCE



Objectives of current research




A. Impact of critical material attributes: crosslinking density




High polymer crosslinking density led to faster drug release



Reduced polymer crystallinity at high crosslinking densities




A. Impact of critical material attributes: polymer functionality




Key determinants?

1. Partition coefficient
2. Diffusion coefficient





A. Impact of critical material attributes: drug physicochemical properties

Evidence of diffusion-controlled release mechanism




Physicochemical attributes dictating drug release from IUSs

Drug	Log P	Mol. Wt (g/mol)	Solubility in polymer (μg/mL)	Melting point (°C)
Medroxyprogesterone	4.1	386.5	7.5730	214.5
Levonorgestrel	3.3	312.4	0.4139	240
Hydrocortisone	1.6	318.1	0.1084	220




A. Impact of critical material attributes: role of additives


Matrix hydrophobicity and drug-excipient interactions are critical



Excipients providing a low-viscosity matrix facilitate faster drug release



Porosity and osmosis may not be significant contributors to drug release




CONTROLLED RELEASE SOCIETY




CRS 2023 ANNUAL MEETING & EXPOSITION
JULY 24-28, 2023 **Paris Hotel » Las Vegas, NV, USA**

THE FUTURE OF DELIVERY SCIENCE


B. Influence of processing parameters



B. Influence of processing parameters



C. Elucidating drug release mechanisms




CONTROLLED RELEASE SOCIETY

CRS 2023 ANNUAL MEETING & EXPOSITION


JULY 24-28, 2023 **Paris Hotel » Las Vegas, NV, USA**

THE FUTURE OF DELIVERY SCIENCE

C. Elucidating drug release mechanisms




Hydrophobic drug in a Hydrophobic Matrix




Hydrophobic drug in a Hydrophilic Matrix


D. Bioequivalence testing through *in vitro* approaches



IUS- 10% w/w drug loading
Saline, 37°C, 100 rpm



IUS- 10% w/w drug loading
45% v/v tert-butanol in PBS, 65°C, 100 rpm



Conclusions

- Tailoring drug release from IUSs □ Controlling polymer crosslinking density
 - Curing conditions
 - Excipient functionality
- Critical attributes □ matrix hydrophobicity, crystallinity, and viscosity.
- Drug release from IUSs □ drug solubility in polymer
 - drug diffusivity in matrix
- Critical processing parameters □ curing temperature and time.
- Accelerated *in vitro* release tests □ quality control and batch testing

Impact

- Propel the development as well as approval of generic equivalent IUSs to advance women's health.
- Tailoring commercial polymers to suit the application.
- Development and optimization of other silicone based controlled release products.



Acknowledgements

Dr. Diane Burgess

Dr. Quanying Bao

Dr. Yuan Zou

Dr. Yan Wang

Kimberly Ma

Kellen Maurus



Burgess Lab members



CONTROLLED RELEASE SOCIETY

CRS 2023 ANNUAL MEETING & EXPOSITION

JULY 24-28, 2023 **Paris Hotel** » **Las Vegas, NV, USA**

THE FUTURE OF DELIVERY SCIENCE