

MERCK
INVENTING FOR LIFE

Tech Session 4: Additive Manufacturing

JUN XU, Ph.D.

Merck: Adam Procopio, Yufang He
Rice University: Jordan Miller Lab, Omid Veisheh Lab

CRS 2022 Annual Meeting & Expo

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada


Advanced Delivery Science




Leveraging Additive Manufacturing to Assay Ultrasound Enhanced Permeability

Goal: To develop a 3D printed *in vitro* vascular model that can be used to assess transport through vasculature for therapeutic applications




Specific Aim 1: To develop a perfused endothelialized vascular model within 3D-printed tissues



Specific Aim 2: To demonstrate this vascular model as an *in vitro* screening tool for assaying ultrasound-induced permeability



Demonstrated a method of consistent endothelialization within 3D printed hydrogels



CRS 2022 Annual Meeting & Expo




Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada






Removal of VEGF resulted in better coverage and more cobblestone morphology

Standard:
Complete Vasculife



Improved:
Vasculife with VEGF
removed




CRS 2022 Annual Meeting & Expo

Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada


Cell-cell junction morphological changes observed within 24 hours of VEGF removal



CRS 2022 Annual Meeting & Expo


Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada




Confirmation of vascular maturity demonstrated by endothelial production of collagen IV

Standard: Complete Vasculife media



GFP HUVEC Collagen IV Hoechst


Improved: Vasculife with VEGF removed



CRS 2022 Annual Meeting & Expo


Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada



Specific Aim 2: To demonstrate this vascular model as an *in vitro* screening tool for assaying ultrasound-induced permeability


Active clinical trials using focused ultrasound



Ying, M. et. al. *Nature Reviews Neurology*. 2021.

Diagnostic (unfocused) ultrasound

- Uses a flat transducer to supply sound waves to a distributed area;
- Relies on low-intensity ultrasound/stable microbubble cavitation to trigger mechanical effects (**sonoporation**);
- Benefit: allows use of lower intensity ultrasound which helps minimize ultrasound bioeffects such as tissue heating;
- **Increasing interest in the use of diagnostic ultrasound for therapeutic applications**



SonoCloud - an implantable, unfocused ultrasound device that provides low intensity ultrasound for blood-brain barrier disruption




CRS 2022 Annual Meeting & Expo

Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada

Sonoporation is utilized to enhance drug delivery through vasculature



Kooiman, K. et. al. *Advanced Drug Delivery Reviews*. 2014.



CRS 2022 Annual Meeting & Expo


Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada




Epifluorescence imaging enabled study of ultrasound-enhanced transport


Experimental set-up



Ultrasound probe positioned on hydrogel



Ultrasound OFF




CRS 2022 Annual Meeting & Expo

Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada

Development of a high-throughput method to quantify vascular permeability changes due to ultrasound

Perfused with 70 kDa Rhodamine-Dextran



fixed region analyzed for fluorescent dye intensity, dI/dt

100 pixels from edge

mean max intensity, I_0

radius, r


automated edge detection

Apparent Permeability Coefficient Calculation

$$P = \frac{dI}{dt} \frac{r}{2I_0}$$

GUI output:
permeability coefficient

	Before US	During US	After US
Whole gel	1.0967	4.5654	5.1727
User ROI	0.8592	2.6132	5.4743
Sliding ROI	1.0930	3.2585	5.7458





CRS 2022 Annual Meeting & Expo

Advanced Delivery Science


July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada

Ultrasound treatment increased apparent permeability coefficient in this hydrogel model




VE-Cadherin-GFP HUVECs provided first real-time acquisitions of endothelial cell-cell junction dynamics under ultrasound

No microbubbles



Microbubbles



CRS 2022 Annual Meeting & Expo

Advanced Delivery Science


July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada




Cell area tracking demonstrated dynamic nature of endothelial cells in response to ultrasound and microbubbles

No microbubbles

Before Ultrasound




After Ultrasound






Microbubbles

Before Ultrasound




After Ultrasound




Conclusion: Developed a new screening platform to help predict therapeutic outcomes *in vitro* before preclinical translation

2D *in vitro* cell culture




3D-printed *in vitro* vascular model

Vascular cross section




3D bioprinted vascular model with HUVECs

Vascular cross section



3D bioprinted BBB model with HBMECs, pericytes, and astrocytes

In vivo animal studies



Adv. Sci. **2020**, 7, 190-3394; *Adv. Therap.* **2021**, 4, 2100154; *Biomaterials Sci.* **2022**, 10, 3158-3173



CRS 2022 Annual Meeting & Expo

Advanced Delivery Science

July 11 – 15, 2022 | Montreal Congress Center, Montreal Canada